




Cxcl11 Antibody

Product Code	CSB-PA09289A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9JHH5
Immunogen	Recombinant Mouse C-X-C motif chemokine 11 protein (22-98AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-X-C motif chemokine 11 (Interferon-inducible T-cell alpha chemoattractant) (I-TAC) (Small-inducible cytokine B11), Cxcl11, Scyb11
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Cxcl11

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA09289A0Rb at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.