



Cxcl11 Antibody, Biotin conjugated

Product Code	CSB-PA09289D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9JHH5
Immunogen	Recombinant Mouse C-X-C motif chemokine 11 protein (22-98AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-X-C motif chemokine 11 (Interferon-inducible T-cell alpha chemoattractant) (I-TAC) (Small-inducible cytokine B11), Cxcl11, Scyb11
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Cxcl11