



TACR1 Antibody, Biotin conjugated

Product Code	CSB-PA09674D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25103
Immunogen	Recombinant Human Substance-P receptor protein (1-29AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Substance-P receptor (SPR) (NK-1 receptor) (NK-1R) (Tachykinin receptor 1), TACR1, NK1R TAC1R
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TACR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.