




GSTP1 Antibody

Product Code	CSB-PA10334A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09211
Immunogen	Recombinant Human Glutathione S-transferase P protein (2-210AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutathione S-transferase P (EC 2.5.1.18) (GST class-pi) (GSTP1-1), GSTP1, FAEES3 GST3
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	GSTP1

Image



Western blot

All lanes: GSTP1 antibody at 2µg/ml

Lane 1: EC109 whole cell lysate

Lane 2: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/15000 dilution

Predicted band size: 24 kDa

Observed band size: 30 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.