



NFKB1 Antibody, Biotin conjugated

Product Code	CSB-PA10354D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19838
Immunogen	Recombinant Human Nuclear factor NF-kappa-B p105 subunit protein (1-961AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear factor NF-kappa-B p105 subunit (DNA-binding factor KBF1) (EBP-1) (Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1) [Cleaved into: Nuclear factor NF-kappa-B p50 subunit], NFKB1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NFKB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.