



GOLM1 Antibody, Biotin conjugated

Product Code	CSB-PA10644D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NBJ4
Immunogen	Recombinant Human Golgi membrane protein 1 protein (36-401AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Golgi membrane protein 1 (Golgi membrane protein GP73) (Golgi phosphoprotein 2), GOLM1, C9orf155 GOLPH2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GOLM1