





MUC16/CA125 Antibody

Product Code	CSB-PA10657A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXI7
Immunogen	Recombinant Human Mucin-16 protein (14067-14328AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:10-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MUC16

Image



Immunohistochemistry of paraffin-embedded human prostate tissue using CSB-PA10657A0Rb at dilution of 1:20



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA10657A0Rb at dilution of 1:20



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.