



S100B Antibody

Product Code	CSB-PA10677A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04271
Immunogen	Recombinant Human Protein S100-B protein (2-92AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein S100-B (S-100 protein beta chain) (S-100 protein subunit beta) (S100 calcium-binding protein B), S100B
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	S100B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.