



S100B Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA10677A0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P04271 |
| Immunogen | Recombinant Human Protein S100-B protein (2-92AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Protein S100-B (S-100 protein beta chain) (S-100 protein subunit beta) (S100 calcium-binding protein B), S100B |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | S100B |