



TDGF1 Antibody

Product Code	CSB-PA10764A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13385
Immunogen	Recombinant Human Teratocarcinoma-derived growth factor 1 protein (32-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Teratocarcinoma-derived growth factor 1 (Cripto-1 growth factor) (CRGF) (Epidermal growth factor-like cripto protein CR1), TDGF1, CRIPTO
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	TDGF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.