



TDGF1 Antibody

Product Code	CSB-PA10769A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13385
Immunogen	Recombinant Human Teratocarcinoma-derived growth factor 1 protein (32-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Teratocarcinoma-derived growth factor 1 (Cripto-1 growth factor) (CRGF) (Epidermal growth factor-like cripto protein CR1), TDGF1, CRIPTO
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	TDGF1