



TDGF1 Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA10769D0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P13385 |
| Immunogen | Recombinant Human Teratocarcinoma-derived growth factor 1 protein (32-150AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Teratocarcinoma-derived growth factor 1 (Cripto-1 growth factor) (CRGF) (Epidermal growth factor-like cripto protein CR1), TDGF1, CRIPTO |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Stem Cells |
| Target Names | TDGF1 |