



COL3A1 Antibody, FITC conjugated

Product Code	CSB-PA11459C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02461
Immunogen	Recombinant Human Collagen alpha-1(III) chain protein (1222-1465AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Collagen alpha-1(III) chain, COL3A1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	COL3A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.