



ALPL Antibody, Biotin conjugated

Product Code	CSB-PA11999D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05186
Immunogen	Recombinant Human Alkaline phosphatase, tissue-nonspecific isozyme protein (28-501AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alkaline phosphatase, tissue-nonspecific isozyme (AP-TNAP) (TNSALP) (EC 3.1.3.1) (Alkaline phosphatase liver/bone/kidney isozyme), ALPL
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ALPL