



PYY Antibody, HRP conjugated

Product Code	CSB-PA12214B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10082
Immunogen	Recombinant Human Peptide YY protein (31-64AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptide YY (PYY) (PYY-I) (Peptide tyrosine tyrosine) [Cleaved into: Peptide YY(3-36) (PYY-II)], PYY
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PYY
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.