



MIB1 Antibody, FITC conjugated

Product Code	CSB-PA13289C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86YT6
Immunogen	Recombinant Human E3 ubiquitin-protein ligase MIB1 protein (5-332AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase MIB1 (EC 2.3.2.27) (DAPK-interacting protein 1) (DIP-1) (Mind bomb homolog 1) (RING-type E3 ubiquitin transferase MIB1) (Zinc finger ZZ type with ankyrin repeat domain protein 2), MIB1, DIP1 KIAA1323 ZZANK2
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	MIB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.