



STEAP1 Antibody, HRP conjugated

Product Code	CSB-PA13354B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHE8
Immunogen	Recombinant Human Metalloreductase STEAP1 protein (3-69AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metalloreductase STEAP1 (EC 1.16.1.-) (Six-transmembrane epithelial antigen of prostate 1), STEAP1, PRSS24 STEAP
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	STEAP1