



gag-pol Antibody, Biotin conjugated

Product Code	CSB-PA15387D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04585
Immunogen	Recombinant Human immunodeficiency virus type 1 group M subtype B Gag-Pol polyprotein protein (1-440AA)
Raised In	Rabbit
Species Reactivity	Human immunodeficiency virus type 1 group M subtype B
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gag-Pol polyprotein (Pr160Gag-Pol) [Cleaved into: Matrix protein p17 (MA); Capsid protein p24 (CA); Spacer peptide 1 (SP1) (p2); Nucleocapsid protein p7 (NC); Transframe peptide (TF); p6-pol (p6*); Protease (EC 3.4.23.16) (PR) (Retropepsin); Reverse transcriptase/ribonuclease H (EC 2.7.7.49) (EC 2.7.7.7) (EC 3.1.26.13) (Exoribonuclease H) (EC 3.1.13.2) (p66 RT); p51 RT; p15; Integrase (IN) (EC 2.7.7.-) (EC 3.1.-.-)], gag-pol
Immunogen Species	Human immunodeficiency virus type 1 group M subtype B
Research Area	Others
Target Names	gag-pol
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.