



rho Antibody, Biotin conjugated

Product Code	CSB-PA16627D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0AG30
Immunogen	Recombinant Escherichia coli Transcription termination factor Rho protein (1-419AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription termination factor Rho (EC 3.6.4.-) (ATP-dependent helicase Rho), rho, nitA psuA rnsC sbaA tsu
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	rho