



rho Antibody, Biotin conjugated

Product Code	CSB-PA16627D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0AG30
Immunogen	Recombinant Escherichia coli Transcription termination factor Rho protein (1-419AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription termination factor Rho (EC 3.6.4.-) (ATP-dependent helicase Rho), rho, nitA psuA rnsC sbaA tsu
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	rho
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.