



Cxcl16 Antibody

Product Code	CSB-PA09337A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8BSU2
Immunogen	Recombinant Mouse C-X-C motif chemokine 16 protein (27-198AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-X-C motif chemokine 16 (Scavenger receptor for phosphatidylserine and oxidized low density lipoprotein) (SR-PSOX) (Small-inducible cytokine B16) (Transmembrane chemokine CXCL16), Cxcl16, Srpsox
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Cxcl16
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.