



Il2rb Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA011650LB01MO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P16297 |
| Immunogen | Recombinant Mouse Interleukin-2 receptor subunit beta protein (24-240AA) |
| Raised In | Rabbit |
| Species Reactivity | Mouse |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Interleukin-2 receptor subunit beta (IL-2 receptor subunit beta) (IL-2R subunit beta) (IL-2RB) (High affinity IL-2 receptor subunit beta) (p70-75) (CD antigen CD122), Il2rb |
| Immunogen Species | Mus musculus (Mouse) |
| Research Area | Immunology |
| Target Names | Il2rb |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |