



Il2rb Antibody

Product Code	CSB-PA011650LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16297
Immunogen	Recombinant Mouse Interleukin-2 receptor subunit beta protein (24-240AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-2 receptor subunit beta (IL-2 receptor subunit beta) (IL-2R subunit beta) (IL-2RB) (High affinity IL-2 receptor subunit beta) (p70-75) (CD antigen CD122), Il2rb
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Il2rb
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.