



DNTTIP1 Antibody, HRP conjugated

Product Code	CSB-PA007089LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H147
Immunogen	Recombinant Human Deoxynucleotidyltransferase terminal-interacting protein 1 protein (1-329AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deoxynucleotidyltransferase terminal-interacting protein 1 (Terminal deoxynucleotidyltransferase-interacting factor 1) (TdIF1) (TdT-interacting factor 1), DNTTIP1, C20orf167 TDIF1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DNTTIP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.