



Thbs2 Antibody, Biotin conjugated

Product Code	CSB-PA023488HD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q03350
Immunogen	Recombinant Mouse Thrombospondin-2 protein (19-232AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thrombospondin-2, Thbs2, Tsp2
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Thbs2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.