



PCLAF Antibody, Biotin conjugated

Product Code	CSB-PA012164DD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15004
Immunogen	Recombinant Human PCNA-associated factor protein (1-111AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PCNA-associated factor (Hepatitis C virus NS5A-transactivated protein 9) (HCV NS5A-transactivated protein 9) (Overexpressed in anaplastic thyroid carcinoma 1) (OEATC-1) (PCNA-associated factor of 15 kDa) (PAF15) (p15PAF) (PCNA-clamp-associated factor), PCLAF, KIAA0101 NS5ATP9 PAF
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PCLAF