



ALB Antibody

Product Code	CSB-PA001561JA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02768
Immunogen	Peptide sequence from Human Serum albumin protein
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum albumin, ALB
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ALB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.