



Lactoferrin Antibody, HRP conjugated

Product Code	CSB-PA00860F0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5MJE8
Immunogen	Sheep Lactoferrin protein (Native Protein)
Raised In	Rabbit
Species Reactivity	Sheep
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lactoferrin
Immunogen Species	Ovis aries (Sheep)
Research Area	Others
Target Names	Lactoferrin