



Lactoferrin Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA00860F0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q5MJE8 |
| Immunogen | Sheep Lactoferrin protein (Native Protein) |
| Raised In | Rabbit |
| Species Reactivity | Sheep |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Lactoferrin |
| Immunogen Species | Ovis aries (Sheep) |
| Research Area | Others |
| Target Names | Lactoferrin |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |