



Lactoferrin Antibody

Product Code	CSB-PA00860E0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5MJE8
Immunogen	Sheep Lactoferrin protein (Native Protein)
Raised In	Rabbit
Species Reactivity	Sheep
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lactoferrin
Immunogen Species	Ovis aries (Sheep)
Research Area	Others
Target Names	Lactoferrin
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.