



ATPAF2 Antibody

Product Code	CSB-PA002425LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N5M1
Immunogen	Recombinant Human ATP synthase mitochondrial F1 complex assembly factor 2 protein (41-289AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP synthase mitochondrial F1 complex assembly factor 2 (ATP12 homolog), ATPAF2, ATP12
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ATPAF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.