



POU2AF1 Antibody, HRP conjugated

Product Code	CSB-PA620973LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16633
Immunogen	Recombinant Human POU domain class 2-associating factor 1 protein (1-256AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	POU domain class 2-associating factor 1 (B-cell-specific coactivator OBF-1) (BOB-1) (OCA-B) (OCT-binding factor 1), POU2AF1, OBF1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	POU2AF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.