



SKAP1 Antibody

Product Code	CSB-PA769799LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86WV1
Immunogen	Recombinant Human Src kinase-associated phosphoprotein 1 protein (1-358AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Src kinase-associated phosphoprotein 1 (Src family-associated phosphoprotein 1) (Src kinase-associated phosphoprotein of 55 kDa) (SKAP-55) (pp55), SKAP1, SCAP1 SKAP55
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SKAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.