



MPO Antibody, Biotin conjugated

Product Code	CSB-PA06899D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05164
Immunogen	Recombinant Human Myeloperoxidase protein (165-278AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myeloperoxidase (MPO) (EC 1.11.2.2) [Cleaved into: Myeloperoxidase; 89 kDa myeloperoxidase; 84 kDa myeloperoxidase; Myeloperoxidase light chain; Myeloperoxidase heavy chain], MPO
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	MPO