





APOC3 Antibody

Product Code	CSB-PA12197A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02656
Immunogen	Recombinant Human Apolipoprotein C-III protein (21-99AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apolipoprotein C-III (Apo-CIII) (ApoC-III) (Apolipoprotein C3), APOC3
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	APOC3


Image



Western blot
All lanes: APOC3 antibody at 2µg/ml + Rat liver tissue
Secondary
Goat polyclonal to rabbit IgG at 1/15000 dilution
Predicted band size: 11 kDa
Observed band size: 11 kDa



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA12197A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA12197A0Rb at dilution of 1:100