



Lgals3bp Antibody

Product Code	CSB-PA012888LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O70513
Immunogen	Recombinant Rat Galectin-3-binding protein (19-574AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Galectin-3-binding protein (Cyp-C-associated protein) (CyCAP) (Lectin galactoside-binding soluble 3-binding protein) (Mac-2-binding protein) (MAC2BP) (Mac-2 BP) (Protein 90K), Lgals3bp
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cancer
Target Names	Lgals3bp
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.