



# Lgals3bp Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA012888LC01RA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O70513
<b>Immunogen</b>	Recombinant Rat Galectin-3-binding protein (19-574AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rat
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Galectin-3-binding protein (Cyp-C-associated protein) (CyCAP) (Lectin galactoside-binding soluble 3-binding protein) (Mac-2-binding protein) (MAC2BP) (Mac-2 BP) (Protein 90K), Lgals3bp
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Research Area</b>	Cancer
<b>Target Names</b>	Lgals3bp