





CNDP2 Antibody

Product Code	CSB-PA856967HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96KP4
Immunogen	Recombinant Human Cytosolic non-specific dipeptidase protein (2-475AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytosolic non-specific dipeptidase (EC 3.4.13.18) (CNDP dipeptidase 2) (Carnosine dipeptidase II) (Epididymis secretory protein Li 13) (Glutamate carboxypeptidase-like protein 1) (Peptidase A), CNDP2, CN2 CPGL HEL-S-13 PEPA
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CNDP2

Image



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA856967HA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA856967HA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.