





# CNDP2 Antibody

<b>Product Code</b>	CSB-PA856967HA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96KP4
<b>Immunogen</b>	Recombinant Human Cytosolic non-specific dipeptidase protein (2-475AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Cytosolic non-specific dipeptidase (EC 3.4.13.18) (CNDP dipeptidase 2) (Carnosine dipeptidase II) (Epididymis secretory protein Li 13) (Glutamate carboxypeptidase-like protein 1) (Peptidase A), CNDP2, CN2 CPGL HEL-S-13 PEPA
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	CNDP2

## Image



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA856967HA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA856967HA01HU at dilution of 1:100