





PKMYT1 Antibody

Product Code	CSB-PA16029A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99640
Immunogen	Recombinant Human Membrane-associated tyrosine- and threonine-specific cdc2-inhibitory kinase protein (900-1101AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane-associated tyrosine- and threonine-specific cdc2-inhibitory kinase (EC 2.7.11.1) (Myt1 kinase), PKMYT1, MYT1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PKMYT1

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA16029A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate tissue using CSB-PA16029A0Rb at dilution of 1:100