



RPA4 Antibody

Product Code	CSB-PA16249A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13156
Immunogen	Recombinant Human Replication protein A 30 kDa subunit protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Replication protein A 30 kDa subunit (RP-A p30) (Replication factor A protein 4) (RF-A protein 4), RPA4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RPA4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.