



CCND1 Antibody, Biotin conjugated

Product Code	CSB-PA004811LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24385
Immunogen	Recombinant Human G1/S-specific cyclin-D1 protein (1-295AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	G1/S-specific cyclin-D1 (B-cell lymphoma 1 protein) (BCL-1) (BCL-1 oncogene) (PRAD1 oncogene), CCND1, BCL1 PRAD1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CCND1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.