



NKTR Antibody

Product Code	CSB-PA015838LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30414
Immunogen	Recombinant Human NK-tumor recognition protein (10-175AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NK-tumor recognition protein (NK-TR protein) (Natural-killer cells cyclophilin-related protein) [Includes: Putative peptidyl-prolyl cis-trans isomerase (PPIase) (EC 5.2.1.8) (Rotamase)], NKTR
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	NKTR