



Transposase for transposon Tn552 Antibody, Biotin conjugated

Product Code	CSB-PA323588LD01FKZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18416
Immunogen	Recombinant Staphylococcus aureus Transposase for transposon Tn552 protein (1-480AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transposase for transposon Tn552 (ORF 480)
Immunogen Species	Staphylococcus aureus
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.