



Thrombin-like enzyme crotalase Antibody

Product Code	CSB-PA520546LA01DYB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	F8S114
Immunogen	Recombinant Crotalus adamanteus Thrombin-like enzyme crotalase protein (25-262AA)
Raised In	Rabbit
Species Reactivity	Crotalus adamanteus
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thrombin-like enzyme crotalase (SVTLE) (EC 3.4.21.74) (Fibrinogen-clotting enzyme) (Snake venom serine protease 2) (SVSP)
Immunogen Species	Crotalus adamanteus (Eastern diamondback rattlesnake)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.