



ASPH Antibody, FITC conjugated

Product Code	CSB-PA002226HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12797
Immunogen	Recombinant Human Aspartyl/asparaginyl beta-hydroxylase protein (75-270AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aspartyl/asparaginyl beta-hydroxylase (EC 1.14.11.16) (Aspartate beta-hydroxylase) (ASP beta-hydroxylase) (Peptide-aspartate beta-dioxygenase), ASPH, BAH
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ASPH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.