



BBE1 Antibody

Product Code	CSB-PA341191LA01EOT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30986
Immunogen	Recombinant Eschscholzia californica Reticuline oxidase protein (24-538AA)
Raised In	Rabbit
Species Reactivity	Eschscholzia californica
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Reticuline oxidase (EC 1.21.3.3) (Berberine bridge-forming enzyme) (BBE) (Tetrahydroprotoberberine synthase), BBE1
Immunogen Species	Eschscholzia californica (California poppy)
Research Area	Others
Target Names	BBE1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.