



EPPK1 Antibody

Product Code	CSB-PA007745LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P58107
Immunogen	Recombinant Human Epiplakin protein (1-225AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Epiplakin (450 kDa epidermal antigen), EPPK1, EPIPL
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	EPPK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.