



LHB Antibody, Biotin conjugated

Product Code	CSB-PA14659D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01229
Immunogen	Recombinant Human Lutropin subunit beta protein (21-141AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lutropin subunit beta (Lutropin beta chain) (Luteinizing hormone subunit beta) (LH-B) (LSH-B) (LSH-beta), LHB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LHB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.