



AHCY Antibody, Biotin conjugated

Product Code	CSB-PA001474HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23526
Immunogen	Recombinant Human Adenosylhomocysteinase protein (2-432AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adenosylhomocysteinase (AdoHcyase) (EC 3.3.1.1) (S-adenosyl-L-homocysteine hydrolase), AHCY, SAHH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AHCY
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.